


# biorettec

## ActivaScrew™ Instrument guide

AO-COMPATIBLE



| Product        | ActivaScrew™<br>Ø 2.0 mm         | ActivaScrew™<br>Ø 2.7 mm         | ActivaScrew™<br>Ø 3.5 mm         | ActivaScrew™<br>Ø 4.5 mm          |
|----------------|----------------------------------|----------------------------------|----------------------------------|-----------------------------------|
| Depth gauge    | B-IS-9020                        | B-IS-9020                        | B-IS-9020                        | B-IS-9021                         |
| Drill bit      | Ø 1.5 mm<br>B-IS-2000            | Ø 2.0 mm<br>B-IS-2700            | Ø 2.5 mm<br>B-IS-3500            | Ø 3.2 mm<br>B-IS-4500 / B-IS-4501 |
| Tap            | Ø 2.0 mm<br>B-IS-2001            | Ø 2.7 mm<br>B-IS-2701            | Ø 3.5 mm<br>B-IS-3501            | Ø 4.5 mm<br>B-IS-4501 / B-IS-4511 |
| Countersink    | B-IS-2004                        | B-IS-3504                        | B-IS-3504                        | B-IS-4504                         |
| Driver type    | HEX shaft, Ø 2.5 mm<br>B-IS-9032 | HEX shaft, Ø 2.5 mm<br>B-IS-9032 | HEX shaft, Ø 2.5 mm<br>B-IS-9032 | HEX shaft, Ø 4.5 mm<br>B-IS-9033  |
| Holding sleeve | B-IS-3503                        | B-IS-3503                        | B-IS-3503                        | B-IS-4503                         |
| Drill Sleeve   | For Ø 2.0 mm<br>B-IS-2005        | For Ø 2.7 mm<br>B-IS-2705        | For Ø 3.5 mm<br>B-IS-3505        | For Ø 4.5 mm<br>B-IS-4505         |
| Large handle   | B-ISC-9013                       | B-ISC-9013                       | B-ISC-9013                       | B-ISC-9013                        |