


bioretec



ActivaScrew™ Cannulated Instrument guide

Product	ActivaScrew™ Cannulated Ø 3.5 mm	ActivaScrew™ Cannulated Ø 4.0 mm	ActivaScrew™ Cannulated Ø 4.5 mm
K-wire	Ø 1.25 mm B-ISC-1250 / B-1251	Ø 1.6 mm B-ISC-1600 / B-ISC-1601	Ø 1.6 mm B-ISC-1600 / B-ISC-1601
Depth gauge	B-ISC-9020	B-ISC-9020	B-ISC-9020
Drill bit	Ø 2.7 mm B-ISC-3500	Ø 3.2 mm B-ISC-4000	Ø 3.5 mm B-ISC-4500
Tap	Ø 3.5 mm B-ISC-3501	Ø 4.0 mm B-ISC-4001	Ø 4.5 mm B-ISC-4501
Countersink	B-ISC-4504	B-ISC-4504	B-ISC-4504
Driver type	HEX shaft for Ø 3.5 mm B-ISC-9032	HEX shaft for Ø 4.0/4.5 mm B-ISC-9033	HEX shaft for Ø 4.0/4.5 mm B-ISC-9033
Holding sleeve	Ø 3.5 mm B-IS-3503	Ø 4.0/4.5 mm B-ISC-4505	Ø 4.0/4.5 mm B-ISC-4505
Drill Sleeve	Double Drill sleeve 2.7/3.5 B-ISC-4505	Double Drill sleeve 4.0/4.5 B-ISC-9034	Double Drill sleeve 4.0/4.5 B-ISC-9034
Large Handle	B-ISC-9013	B-ISC-9013	B-ISC-9013